

SARS-CoV-2

Proficiency Testing Scheme



NEXT CALL

Due to the necessity to evaluate the evolution of the global pandemic, the development of detection test is being prioritized not only in patients but also in the environment.



In this regard, is very important to demonstrate the reliability of the results obtained, and the PT Schemes are a fundamental tool to evaluate the technical competence of the assays developed.

ielab will open very soon the registration for a **PT Scheme** about the **SARS-CoV-2 (COVID-19)** detection by means **RT-qPCR**.

SAMPLES

Samples employed will be synthetic, no-infectious, and will content RNA from several SARS-CoV-2 target genes (E, N, ORF1ab, RdRP y S), which will allow evaluating the complete virus detection methodology. The fields of application are:

- ❖ **CLINICAL / SANITARY**
- ❖ **ENVIRONMENTAL:** Waste water, Reclaimed water, Air quality, etc.
- ❖ **SURFACES**



RESULTS

Results will be reported qualitatively (Presence/Absence).

EXPERIENCE

ielab is a company belonging to SUEZ group that has a wide experience providing services and products related to the application of quality in laboratories. **ielab** is **accredited** by ENAC as a PT Scheme provider under the **ISO 17043** standard, with more than 15 years of experience. In the international PT Schemes organized by **ielab**, more than 1450 laboratories distributed all over the world participate annually.

If you need more information please, do not hesitate to contact us at the e-mail address comercial@ielab.es, we will be glad to assist you.

ielab

T +34 966 10 55 01

C/ Dracma 7 - P. I. Las Atalayas
03114 Alicante (España)

www.ielab.es